

Fragment Length Analysis

Sample Amounts and Concentrations for Plates

Volume	Pipetting Scheme
Each sample should have a volume of at least 10 µl	Arrayed in a 96-well plate

Sample Preparation

Labeled PCR Products

- Dilute your PCR product with Hi-Di formamide to reach the desired signal intensity (e.g. 1:2 - 1:100).
- Pipette 10 µl of the diluted, labeled PCR products into a 96-well plate.

Ready-to-load

- Dilute your size standard 1:100 with Hi-Di formamide
- Distribute 10 µl size standard/formamide mix into a 96-well plate (PCR plate, semi-skirted, ABI style).
- Dilute your PCR product to reach the desired signal intensity (e.g. 1:2 - 1:100).
- Add 1 µL of diluted PCR product to 10 µL size standard/formamide mix.

Remark: Safe shipment of liquid solutions requires good sealing of your 96-well plate. It is recommended that you seal your 96-well plate with 8-cap stripes, 96-cap mats or heat sealing.

Matrix, Filter Sets and Dye Combinations

Dye Set Name	Matrix Type	Blue	Green	Yellow	Red	Orange
Powerplex_4C	4 dye	FAM	JOE	TAMRA	CXR*	
Powerplex_4C-RCT	4 dye	FAM	JOE	TAMRA	CXR*	
DS-30_D	4-dye	6-FAM	HEX	NED	ROX*	
DS-30_D-RCT	4-dye	6-FAM	HEX	NED	ROX*	
DS-02_E5	5 dye	dR110	dR6G	dTAMRA	dROX	LIZ*
DS-02_E5-RCT	5 dye	dR110	dR6G	dTAMRA	dROX	LIZ*
DS-33_G5	5 dye	6-FAM	VIC	NED	PET	LIZ*
DS-33_G5-RCT	5 dye	6-FAM	VIC	NED	PET	LIZ*
Mic-G5	5 dye	FAM	ATTO532	ATTO550	ATTO565	Dyomics630*
Mic-G5-RCT	5 dye	FAM	ATTO532	ATTO550	ATTO565	Dyomics630*

RCT: This dye set reduces signal intensity and potential crosstalk of capillaries. Higher concentration peaks can be used without going offscale.

** Dyes are reserved exclusively for the size standard and cannot be used for labeling fragments.*

Order Form Completion for Plates

Prior to shipping your samples to Microsynth, please proceed as follows to complete your order:

1. Enter our webshop on www.microsynth.com (click on "LOGIN SHOP")
2. Click on „**Genetic Analysis**“ in the green Analysis Services area
3. Click on one of the options under **Fragment Length Analysis**.
4. Fill in the order form and submit your order
5. Pack your samples into a transparent plastic bag (important: one bag per order) **with** the printed **order form**.
6. Drop your sample package into the nearest Microsynth sample drop box (if available in your vicinity) or alternatively use mail shipment.

Need More Information?

Microsynth AG

Schützenstrasse 15
9436 Balgach
Switzerland
Phone: +41 71 726 10 04
Email: sanger.support@microsynth.ch

Microsynth Seqlab GmbH

Maschmühlenweg 36
37081 Göttingen
Germany
Phone: +49 551 37 000 15/17
Email: info@microsynth.seqlab.de

Microsynth Austria GmbH

Leberstrasse 20
1110 Wien
Austria
Phone: +41 71 726 10 04
Email: sanger.support@microsynth.ch